Weldix® Fadenanker
Weldix® Suture Anchor

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Contents

1	General Information	15
2	Intended Use	15
2.1	Term definition.....	15
2.2	Compatibility.....	15
2.3	Notes	16
3	Warranty.....	16
4	User's Inspection	16
5	Hotline	17
6	Users.....	17
7	Implant Materials	17
8	Indications	17
9	Contraindications.....	18
10	Potential Side Effects and Complications.....	19
11	Packing of the Implants.....	19
12	Structure and function of the implant system	20
12.1	Weldix® Anchor	20
12.2	Twist Drills.....	20
13	Storage	21
14	Intraoperative Procedure	21
15	Postoperative Procedure.....	21
16	Disposal	22

1 General Information

These instructions for use apply to the product line of:

VetWelding AG
Muehlebach 2
CH-6362 Stansstad/NW - Switzerland
Tel. +41 41 530 70 99

info@vetwelding.com – <https://vetwelding.com/>

Non-compliance with the instructions for use can lead to injuries of the patient!

Please note:

- Every user is required to read these instructions for use completely and follow them, and must inform the patient about the usage instructions, the safety warnings and the safety instructions contained in these instructions for use.
- Keep the Instructions accessible to all users at all times.

2 Intended Use

The Weldix® Anchor is intended for sutures or tissue fixation in the foot, ankle, knee, wrist, elbow and shoulder. The Weldix® anchor is exclusively intended for insertion with the BoneWelder Vet.

2.1 Term definition

The Weldix® anchor stands for implants whose base material consists of biocompatible and fully bioresorbable polylactide (PLDLLA 70:30).

2.2 Compatibility

The Weldix® anchor may only be used with the thread types, materials and sizes intended for this purpose.

Suture type	Suture size
Monofilament, resorbable	USP #2-0 to #2
Polyglyconate (copolymer of glycolic acid /trimethylene carbonate copolymer)	
Monofilament, non-resorbable	USP #2-0 to #1
Polypropylene Polyamide 6 and 6.6	
Multifilament, non-resorbable	USP #3-0 to #2
Polyester (polyethylene terephthalate)	
UHMWPE/polyester	USP #4-0 to #2

2.3 Notes

- In accordance to their characteristic attribute as resorbable implant system, Weldix® implants are bio-degradable and excreted from the patient's body. Therefore, the implants cannot be reused.
- Reprocessing is not possible because the product properties would be damaged.
- The VetWelding implant systems are intended for surgical procedures in which an internal fixation by biologically degradable implants is required for aligning, reconstructing and stabilizing bone tissue.
- Different surgical methods and materials can be used for the application of the implants. Type, structure and application of each implant system are determined by the user after evaluation of the patient's anatomical and functional requirements and are therefore the responsibility of the user (treating Veterinarian or surgeon).
- The product has to be used according to indication. VetWelding assumes no responsibility for complications caused by incorrect indication, incorrect choice of implant, incorrect combination of system components and surgical technique as well as limitations of the treatment method or insufficient hygiene.
- All implant components, instruments and equipment must be in sterile, ready-to-use condition before they are used.
- All absorbable implants are sterilized in the factory and packed in sterile packing.
- The surgeon must be familiar with the Weldix® material as well as with the corresponding system components BoneWelder® Vet. It is assumed that their specific instructions for use are known. The user is familiar with the surgical procedure and is responsible for the selection of the correct implants and their correct insertion.

WARNING

Combination of different materials or systems, with the exception of compatible sutures, is not permitted!

VetWelding implant systems must not be combined with systems of other manufacturers.

Do not apply instruments which are not specifically determined for Weldix® and BoneWelding® technology. Doing this can lead to failure of the implant system.

3 Warranty

Our Standard Terms and Conditions of Sale effective at the time shall apply. Agreements diverging from these Standard Terms and Conditions do not restrict the legal rights of the buyer.

4 User's Inspection

Immediately upon receipt, the goods must be checked for completeness and potential damage in transit. Potential transport damage and knocking over the thermo indicator can have negative influence on the quality and the working condition of the specific implant.

Damaged products must be disposed of. They must not be used in surgical procedures

5 Hotline

Should you have any questions on how to handle the unit/product or question on its clinical application, please do not hesitate to contact us.

Tel: + 41(0) 41 530 70 99 (CET / GMT+1 8:30 – 17:30)

Tel: + 41(0) 76 331 65 65 (CET / GMT+1 8:00 – 24:00)

NOTICE

Each package comes with a specific lot number. Please always indicate this number together with the item number in case of complaint.

6 Users

- The implant system from VetWelding may only be used by competent medical staff who, due to their specialized training, experience and knowledge of the relevant regulations have the necessary qualification for using the system correctly.
- Users must have adequate knowledge of the rules of medical practice, as well as the state of science and the contents of pertinent scientific articles published by medical authors.
- It is the user's responsibility to avoid or minimize the general risks posed by surgical interventions.
- The user
 - bears the responsibility for the proper performance of the surgical procedure;
 - must be thoroughly familiar with the surgical technique as well as with the set of implants and the instruments and accessories used;
 - is responsible for selecting the right implant components and implanting them properly.

7 Implant Materials

The Weldix® Anchor -Polymer is a 100% poly-D,L-lactic acid (PDLLA). PDLLA is well established as an absorbable polymer for use in veterinary medicine.

The core of the degradation process: The complex polymer chains absorb the water contents (H₂O molecules) of surrounding body fluids through a process called "hydrolysis". The stored water initiates the degradation process by continuously breaking down the long polymer chains into ever shorter structures or simpler molecules. Metabolic pathways subsequently transform the molecules into carbon dioxide and water; both of these compounds are discharged naturally.

8 Indications

The VetWelding's Weldix® implant system is intended for surgical procedures where internal fixation by resorbable implants is required for reconstruction and stabilization of bone tissue (joints), such as:

- Hip dislocation (capsular repair, repair of the ligamentum teres femoris).
- Repair of the tarsal collateral ligament

- Achilles tendon repair
- CCL extracapsular repair
- Knee collateral ligament repair
- Elbow/shoulder dislocation repair

For more details see the respective product documentation.

9 Contraindications

Non-observance of the contraindications poses a danger of serious or even lethal patient injury!

The veterinarian surgeon's education, training and professional judgment must be relied upon to choose the most appropriate device and treatment.

Conditions presenting an increased risk of failure include:

- Implantation in bone areas with minimal or no cancellous bone structure. The BoneWelding® procedure requires cancellous bone for implant integration. The device must not be implanted in the absence of cancellous bone.
- Active or latent soft tissue and/or skin infections at the site of the operation.
- Local inflammation of any origin at the repair site.
- Hypersensitivity to foreign bodies.
- Suspected material sensitivity or allergies and/or reactions to the implant materials used.
- Autoimmune diseases.
- Inadequate circulation or other pathological changes of the bone and/or soft tissue at the repair site.
- A bad general state of health and/or vascular disorders and/or metabolic disorders (such as diabetes) that may affect healing.
- Bone tumors in the area of the implant anchoring.
- Insufficient quantitative or qualitative hard and soft tissue.
- Concomitant diseases, e. g. degenerative disease processes with negative influence on successful healing.
- Osteoporosis or osteomalacia and other serious damage to the bone structure conflicting with solid anchoring of the implant components.
- Factors which may also impair the success of operations include
 - Incorrect implantation technique.
 - Use of a non-compatible suture material (e.g. PDS).
 - Inadequate (too little, too short) postoperative load reduction and immobilization for healing.

10 Potential Side Effects and Complications

- Implanting foreign bodies can generally cause inflammatory or allergic reactions.
- Decrease of bone density and/or bone necroses.
- Implant loosening, tear-out or fracture.
- Vascular changes.
- Luxations and subluxations.
- Early and late infections, surface and/or deep infections.
- Damage to nerves and vessels.
- Hematomas and wound healing dysfunctions.
- Periarticular calcification.
- Restricted joint function and movability.
- Restricted joint strain and joint pains.
- Allergic reaction or sensitivity to the implant material.
- Nerve damage due to surgical trauma.
- Insufficient bone formation or delayed bone healing, osteo- porosis, osteolysis, osteomyelitis, slowed re-vascularization or infections, bone necrosis that can cause a loosening, bending, tearing or breaking of the implant systems.
- Pain, discomfort, abnormal sensations or touch ability due to the implants.
- Increased reaction of the soft tissue in proximity to the osteotomy and/or the implant systems.
- Insufficient healing.
- During the absorption phase a temporary increase of volume (swelling) of the implants can occur.

Besides the above-mentioned side effects, complications such as infections, nerve damage and pains that are not connected to the implant system can occur in any surgical procedure.

11 Packing of the Implants

- Resorbable implants from VetWelding will generally be delivered in sterile conditions (gamma ray sterilization). The implants with the label **STERILE R** can be taken out of their packing and applied immediately without further pre- treatment.
- After the delivery of the product, the sealing of the original packaging needs to be checked.
- The temperature indicator must immediately be checked on receipt. Products whose thermo indicator has turned from silver to black needs to be rejected. It must not be used.
- The expiration date must also be checked. Symbol: . It is indicated on the outer protective packaging.
- The transport/storage packing of sterile implants is unsterile. This is why sterile implants must be unpacked from their transport packing outside the sterile area of an operation room and given into the sterile area without compromising the product's sterility.
- Re-sterilization of the implants is not allowed.

- Implants coming in damaged or opened sterile packages and implants with passed expiry date have to be considered as being unsterile. Therefore, they have to be rejected.

NOTICE

Each package comes with a specific lot number. Please always indicate this number together with the item number in case of complaint.

The implant components used can be documented in the patient file with their item numbers, item descriptions and batch and/or serial numbers. Only the documentation of the serial numbers ensures a clear traceability of the implants used.

12 Structure and function of the implant system

12.1 Weldix® Anchor

CAUTION

System failure due to deformation of Weldix® Anchor!

Weldix® Anchors must not be heated or bent before they are inserted correctly.

- Weldix® Anchors are only approved for use in combination with the BoneWelder® Vet system (REF 01-00-001). For detailed information on the BoneWelder® Vet read the Instructions for Use that are supplied with the BoneWelder® Vet.
- When using the Weldix® system, make sure that the core hole of the Weldix® Anchor is only pre-drilled with special drills. Tapping is not required.

12.2 Twist Drills

- When twist drills are used, proper cooling must be provided in order to avoid heat damage and necroses on the bone.
- Follow the manufacturer's instructions regarding the handpiece that is used with the twist drill. The maximum number of turns per minute has to be respected.
- Too extreme application of force can lead to excessive strain that can damage or break the drill.
- Take protective measures so that a potential breaking of the twist drill does not injure the patient, the user or others.
- In order to avoid possible damage to tissue close to the drill, one has to work with respective drill protection (e. g. drill cannula).
- The drills must be used with a maximum speed of 1,000 rpm. The drill must be positioned orthogonally. Start drilling using medium turns per minute then increase the number of turns applied.
- If you notice that the drill is imbalanced, the drilling process must be stopped immediately. The defective drill must be replaced.

13 Storage

The implants must be stored at room temperature.

A maximum temperature of 98.6°F (37°C) must not be exceeded.

- Do not use the implant material any more, if the thermo indicator has changed from silver to black.

14 Intraoperative Procedure

The product may only be implanted with the specially designed drill bit (REF 01-05-001 and REF 01-05-002) and the specially designed sonotrode (REF 01-04-001). Weldix® implants may only be implanted with the Beone-Welder® Vet Ultrasound System (REF 01-00-001) provided for this purpose.

- The implant should be pushed onto the sonotrode as far as possible.
- Load with the desired and suitable suture. Position the suture into the groove at the anchor's distal end.
- One or two threads can be used per anchor.
- Pre-drill the bone by using the drill bit (REF 01-05-001 und REF 01-05-002).
- Drill until the drill stop reaches the bone surface.
- Push the anchor into the countersink of the pre-drilled hole before activation the ultra-sound.
- Keep the anchor pointed straight into the pre-drilled hole during the entire insertion.
- No tension on the suture during the insertion.
- Insert the anchor in one continuous movement.
- Insert the anchor in less than 3 seconds.
- Keep pressing the anchor into the drill hole as you activate the ultrasound (pushing the activation button on the handpiece) until the anchor sits flush with the surface.
- Stop the ultrasonic process (release activation button on the handpiece) as soon as the implant is fully inserted into the bone to avoid excessive heat development
- Allow the implant material to completely solidify for five seconds before removing the handpiece.
- After the Weldix Anchor® is fixed in the bone, gently pull the suture strands to confirm adequate bone fixation and no suture damage.

15 Postoperative Procedure

In the resorption phase, a temporary increase of volume (swelling) of the implants can occur. This must also be considered intraoperatively.

⚠ WARNING**Danger of injury and risk of system failure due to improper behavior of the patient!**

Inadequate postoperative load reduction and immobilization for healing can lead to system failure or injuries of the patient.

Therefore, please note:

- Instruct owner and give postoperative care.
 - All physical activities that can affect the implant system or the bone elements that are connected to it must be refrained during the entire duration of the treatment. This applies especially to activities that involve impacts or blows.
 - Suitable measures must be taken to prevent overloading of the anchorage.
-

16 Disposal

The national disposal regulations and guidelines must be observed when disposing packaging material and potentially infectious material (e. g. after removing metal).

VetWelding AG

Muehlebach 2

CH-6362 Stansstad/NW- Schweiz

Tel. +41 41 530 70 90

info@vetwelding.com – <https://www.vetwelding.com>

A Company of WoodWeldingSA using BoneWelding@Technology

The BoneWelding@technology is covered by the intellectual property rights of and licensed from WoodWeldingSA, Switzerland

For trademark information please refer to www.vetwelding.com/about-us

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Copyright by VetWelding AG · Alle Rechte vorbehalten ·